Manugrip

Dear Customer, We appreciate the confidence placed on a product EMO. Choosing a brace Manugrip brand you have purchased a quality product and high level medical doctor. Read these instructions carefully. If you have any questions please contact your doctor or your specialist.

REGULATIONS

The medical device complies with Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices and the respective national regulations. Class I medical device, non-sterile, without measurement function.

DESCRIPTION AND FEATURES

The MQ701 Manugrip is an immobilizing brace wrist joint. It has an anatomical shape that stabilizes and immobilizes the joint maintaining the mobility of the fingers. Technologically Advanced: manufacture by micro - fusion avoids overpressure and scratches as it is made without seams. Lightweight : materials made of maximum strength and minimum thickness. It incorporates aluminum plates and removable shaped: palmar, ulnar and radial . Placement with a single hand: adjustment straps close on themselves to facilitate colocation.El adjustment is made using some straps passing through buckles mutiplican a pulling force and easy adjustment. Breathable fabric foamed. Bilateral. S (14-18 cm), L (18-22 cm)

INDICATIONS

Injuries to the wrist -joint without fracture, controlled immobilization of the joint, joint instability, carpal tunnel syndrome, post -operative synovitis,degenerative osteoarthritis, recurrent sprains wrist.

INSTRUCTIONS FOR USE AND PLACEMENT

The product should be adapted by the prosthetist or by the physician individually for each patient. The correct choice of size is essential to achieve maximum effectiveness of treatment and prolong the life of the garment. Loosen the straps and secure about themselves . Check the conformation of the strips and adapted to the morphology of the patient , if necessary. Enter the injured hand and adjust. The patient can get on and off the cuff without mobilizing Just The Tip, facilitating rehabilitation exercises and personal hygiene.

PRECAUTIONS

Follow the manufacturer's directions for proper orthosis placement. Follow the instructions or recommendations of the prescriber or healthcare personnel. For a personalized adaptation, consult

a specialist. If there is ongoing pain, it is advisable to contact a healthcare professional. If you experience any allergic reaction to remove the product and consult with the prescribing physician. Do not expose the garment to a flame or source of heat at very high temperatures. Store at room temperature, preferably in the original container. Warranty within the framework of the laws.

CLEANING AND MAINTENANCE

Hand wash with warm water (max. 30 °) and neutral soap. During washing, if the garment has Velcro, they must stick together or keep them closed. Dry at room temperature. Do not iron the garment shop. Contact with ointments, creams, alcohol and solvents can damage the garment.


QUALITY CERTIFICATES

All materials used have been tested , comply with European regulations, are high quality and offer maximum safety and comfort during use. Öko -Tex Standard 100: textile eco-label awarded after the completion of a catalog of criteria for analysis that guarantees the absence of substances harmful to health. Certifies the absence of harmful substances in all materials with which a garment is made.

NOTICE

To dispose of the product, use the solid waste container. For the disposal of the container, depending on its material, use the paper / cardboard container or the plastic container. In any case, you must strictly comply with the regulations of the country. Any serious incident related to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user and / or the patient are established.

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